

## Shanghai, China ETO

### Facility Address:

333, Shuang Hui Road  
201 308 YangShan Free Port, Shanghai, P.R.. China

### Facility Services Capabilities:

Ethylene Oxide sterilization  
EOStat® Rapid Processing  
Laboratory – B.I. Testing, Product Bioburden  
Cycle Development & Process Validation,  
Logistics



### Process Capability:

4 Sterilizers Maximum pallet height 2,4m Pallet size 1.2m x 1.0m or 1.2m x 0.8m

### Key Contacts for:

Sr. Director  
Asia Pacific Plant Manager  
Customer Service  
Logistics  
Sales

### Contact Name:

Erik Jia  
Jessica Yuan  
Amei Xu  
Blake Tan  
Clare Hu

### Telephone:

+86 21 6828 0017  
+86 21 6828 0157  
+86 21 6828 0237  
+86 21 6828 0216  
+86 21 6828 0389

### E-mail:

ejia@sterigenics.com  
jyuan@sterigenics.com  
axu@sterigenics.com  
btan@sterigenics.com  
chu@sterigenics.com

### Hours of Operation:

Weekdays  
Saturday / Sunday  
Holidays Closure Dates

### Processing:

24 hours  
24 hours

### Delivery/Collection:

9:00am–9:00pm  
9:00am–9:00pm

### Office Inquiries:

9:00am–5:00pm  
Closed  
National Holidays

### Local Hotels:

Howard Johnson Hotel  
Crowne Plaza

+ 86 (21) 3828 8888  
+ 86 (21) 2033 9999 x8622